

User Manual



ImplantCenter M+



Contents

1 Documentation	7
1.1 Associated documentation	7
1.2 Electronic documentation	7
2 Required information	9
2.1 Indication for use	9
2.2 Operating principle	9
2.3 Date of first inclusion of EC marking	9
2.4 Latest document update	9
2.5 Repairing or modifying the medical device	9
2.6 Warranty	10
2.7 Accessory usage conditions	10
3 Removal from packaging, installation, connections	11
3.1 Removing the medical device from its packaging	11
3.2 Installing the medical device	11
3.3 Installing cords	11
3.4 Connecting the medical device to the electrical network	11
3.5 Fixing the medical device to a non-removable support	12
4 Medical device description	13
4.1 Control unit	13
4.2 Handpieces and I-Surge LED micromotor	13
4.3 Cords	13
4.4 Inlets	13
4.5 Switch	13
4.6 Mains connector	13
4.7 Fuse recess	13
4.8 Pump	13
5 User interface	14
5.1 I-surge mode	14
5.2 Piezotome® M+ mode	15
5.3 Newtron® mode	16
6 Medical device use	17
6.1 Preparing the medical device	17
6.1.1 Connecting the medical device to the electrical network	17
6.1.2 Footswitch	17
6.1.3 Installing the footswitch	17
6.1.4 Connecting the I-Surge LED micromotor	17
6.1.5 Connecting the Piezotome® M+ handpiece	17
6.1.6 Connecting the Newtron® handpiece	17
6.1.7 Installing the irrigation system	17
6.1.8 Purging the medical device	18

6.2 Switching on the medical device	18
6.3 Selecting the operating mode	18
6.3.1 Common operations using the footswitch	18
6.4 Connecting and disconnecting M+ accessories during use	19
6.5 Switching off the medical device	19
7 Configuring the medical device	21
7.1 Adjusting the medical device volume	21
7.2 Adjusting the screen brightness	21
7.3 Adjusting the lighting timer	21
7.4 Recovering the medical device's factory settings	22
7.4.1 I-Surge mode factory settings	22
7.4.2 Piezotome mode factory settings	22
7.4.3 Newtron mode factory settings	22
7.5 Personalising a contra-angle or handpiece	22
8 Cleaning, disinfecting and sterilising	23
8.1 Cleaning and disinfecting the medical device	23
8.2 Cleaning, disinfecting and sterilising accessories	23
9 Monitoring and general maintenance of the medical device	25
10 Maintenance	27
10.1 Touch-sensitive screen messages	27
10.1.1 Problem at start up	27
10.1.2 Missing handpiece	27
10.1.3 I-Surge LED micromotor missing	27
10.1.4 Finding the software version	28
10.2 Identifying incorrect operation	28
10.2.1 Not working	28
10.2.2 The power is not as expected	28
10.2.3 No spray or very little amount of spray	29
10.2.4 Light not working or too dim	29
10.2.5 Water leakage	29
10.2.6 Ultrasounds not working	29
10.3 Corrective maintenance	30
10.3.1 Replacing the fuses	30
11 Technical specifications of the medical device	31
11.1 Identification	31
11.2 Control unit	31
11.3 Generator	31
11.4 Screen	31
11.5 Length of cords	31
11.6 I-Surge LED micromotor	32
11.7 Irrigation	32
11.8 Footswitch	32
11.9 Environmental characteristics	32
11.10 Environmental restrictions	32
11.11 Main performance characteristics	32
12 Regulations and standards	33
12.1 Official texts	33
12.2 Medical class of the device	33
12.3 Symbols	33
12.4 Manufacturer identification	34
12.5 Branch addresses	35
13 Disposal and recycling	37
14 Index	39

15 Glossary	41
-------------------	----

Foreword

The SATELEC[®] medical device that you are about use is designed for professional use only. It is therefore a key tool with which you will provide treatment within the context of your work.

These medical devices are designed to be used exclusively within a hospital or private clinic operating theatre. Patients will be under general anaesthesia unless the procedure only requires local anaesthesia.

To ensure optimum safety for yourself and your patients, comfort in your daily practice and to benefit fully from your medical device's technology, please read the documentation provided carefully.

If you have received this medical device by mistake, please contact the supplier to arrange for it to be collected.

Please refer to the instructions for the entire range of SATELEC[®] M+ ultrasonic medical generators for information about the following:

- documentation format;
- documentation archiving period;
- warnings concerning user and patient populations;
- treatment area;
- medical device usage interactions, contraindications and prohibitions;
- electromagnetic compatibility;
- disposal and recycling of the medical device;
- manufacturer responsibility.

Please refer to the various cleaning, disinfection and sterilisation protocols for information about the following:

- preparation of parts for sterilisation;
- detailed manual and automatic protocols;
- information concerning conditioning for sterilisation;
- recommendations for the inspection of parts.

1 Documentation

This document contains the following information:

- indications for use
- medical device description
- installation of the medical device
- medical device use
- preparation prior to cleaning and disinfecting the medical device
- monitoring and general maintenance of the medical device
- maintenance to be performed by the user.

1.1 Associated documentation

This document must be used in association with the following documents:

Document title	References
General instructions relating to the complete range of M+ medical ultrasonic generators	J57821
General instructions relating to all M+ medical handpieces	J12861
Method for consulting electronic user instructions	J00000
Cleaning, disinfection and sterilisation protocols for the Piezotome [®] M+ handpiece	J12801
Cleaning, disinfection and sterilisation protocols for the maintenance-free I-Surge LED implantology motor	J28721
I-Surge LED micromotor Quick Clean	J27120
User manual for I-Surge LED micromotor	I27210
ImplantCenter [™] M+ Quick Start	J27260
ImplantCenter [™] M+ Quick Clean	J27261
ImplantCenter [™] M+ user manual	J27251
Piezotome [®] M+ handpiece user manual	J57521
Newtron [®] handpiece user manual	I12610
Newtron [®] handpiece Quick Clean	J12930
Pack M+ Quick Start	J57870
Pack M+ Quick Clean	J57871
Pack M+ user manual	J57841

1.2 Electronic documentation

The user instructions for your device are provided in electronic format and not in printed format. However, you can request a free printed copy of the user instructions within 7 days via our website, by telephone or in writing.

The electronic user instructions are available in PDF format (Portable Document Format) and you will need to have a PDF file read software installed to read the instructions.

The device user instructions can be consulted at the following address:

www.satelec.com/documents



It is important for you to have read and understood the content of the user instructions relating to the use of your device and its accessories prior to use.

We recommend that you visit the website regularly to consult and/or to download the latest version of your device's user instructions.

2 Required information

2.1 Indication for use

ImplantCenter™ M+ is a surgical device that combines an M+ surgical or ultrasonic piezoelectric handpiece and a rotary surgical motor.

The Piezotome® M+ piezo-ultrasonic surgical handpiece can be used in combination with M+ tips for cutting into bone, metal and bone substitutes.

The device's Piezotome® M+ surgical piezoelectric part can be used for osteotomy, osteoplasty, decortication, drilling, shaping and smoothing of teeth and bone in various surgical procedures, including but not limited to general, orthopaedic, otolaryngological, maxillofacial, oral, hand and foot, neurosurgery, spinal and plastic/reconstructive surgery.

The Newtron® piezo-ultrasonic handpiece associated with dental tips can be used for mechanised ultrasonic dental treatments, e.g. by periodontics, endodontics, scaling and prosthetics practices.

The rotary surgical I-Surge LED motor is indicated for use in dentistry such as implantology, dental surgery and endodontics as well as maxillofacial surgery.

The device can be used on patients of all ages, including paediatric patients.

2.2 Operating principle

An electrical signal emitted by the medical device is supplied to the dental piezo-ultrasonic handpiece. It comprises a piezoelectric ceramic transducer, which converts the electrical signal into ultrasonic vibrations. Mechanical vibrations are transmitted to a tip attached to the end of the ultrasonic handpiece.

The vibrations, applied with a very gentle pressure, cut the bone with ultrasounds. They act distinctively on the bone, minimising the risks of damaging the soft tissue, thereby guaranteeing a more precise and safe procedure with less strain on the surgeon's hand. This ultrasonic piezoelectric surgical procedure also improves bone healing.

The I-Surge LED rotary surgical micromotor can be used in combination with a handpiece, a rotary contra-angle and other rotary drills or oscillating saws. It can be used for cutting into bone, placing screws, broaches and retaining wires.

The type of coupling of the I-SURGE LED micromotor developed by SATELEC® is compatible with most contra-angles and straight handpieces.

| Contra-angles and straight handpieces are not supplied by SATELEC®

2.3 Date of first inclusion of EC marking

2013

2.4 Latest document update

01/2016

2.5 Repairing or modifying the medical device

Contact the supplier of your device. Using the services of an unapproved repairer could render your device dangerous for you and your patients.

Do not repair or modify the device without seeking the prior permission of SATELEC®.

If the device is modified or repaired, specific checks and tests must be carried out to ensure that the medical device is still safe to use.


In the event of doubt, contact an approved dealer or the SATELEC® customer service team:

www.acteongroup.com

satelec@acteongroup.com

SATELEC® at the request of technical personnel working for the network of dealers approved by SATELEC®, provide all information required to repair the faulty parts on which they may perform repairs.

2.6 Warranty



The user may not remove any of the screws shown on this view, as this would invalidate the medical device's warranty.

2.7 Accessory usage conditions

Accessories, handpieces and motors must be cleaned, disinfected and sterilised prior to use.

3 Removal from packaging, installation, connections

3.1 Removing the medical device from its packaging

When you receive your medical device, check for any damage that may have occurred during transportation. If you have received this medical device by mistake, please contact the supplier to arrange for it to be collected.

If you have any questions or requirements, contact your supplier.

The ImplantCenter™ M+ includes the following items:

- an ImplantCenter™ M+ unit
- an I-Surge LED implantology motor and cord
- a multifunction footswitch
- two peristaltic pump recesses
- two Irrigation solution brackets
- five irrigation lines with 10 single-use sterile clips
- two handpiece supports
- 30 single-use, sterile irrigation line perforators
- a power cord
- a [J27260] Quick Start guide, a [J27261] Quick Clean guide
- an attachment kit.

3.2 Installing the medical device

Place the control unit in the position that is suitable for your activity. Check that the cords do not hinder the movement or free circulation of anyone. The medical device must be placed on a secure and flat surface or a surface with a maximum slope of 5 degrees.

Fix your medical device using the attachments provided to ensure that the device cannot be removed without the use of a tool.

Adjust the position of your medical device to correspond to your angle of vision and the characteristics of your workstation, e.g. lighting or distance between the user and the medical device.

Ensure that your medical device is readily accessible.

| Do not install your medical device near or on another device.

3.3 Installing cords

Check that the cords do not hinder the movement or free circulation of anyone. Make sure that it is not possible to wheel over or walk on the different cords. Cords attached to their handpiece or I-Surge LED micromotor must be readily accessible. Make sure that the cord is slack during use.

Never rotate the handpiece connector on its cord as this can damage your medical device.

| Never wrap the handpiece cord around the medical device. Do not put the medical device cords in a cable cover or a cable tray.

3.4 Connecting the medical device to the electrical network

Switch the medical device OFF (position 0) and check that the mains voltage is compatible with that indicated on the medical device. Next, connect the cord to the wall socket in compliance with the standards in force in the country of use.

The medical device is equipped with a protective ground connection and must be connected to a mains power supply with a protective ground.

Do not plug the medical device into an extension lead and do not put the mains cord in a cable cover or cable tray.

A different voltage would cause damage to the medical device and could injure the patient and/or user. Any variation in the electrical network voltage or electromagnetic field that is non-compliant with the limits in force, could interfere with the medical device's operation.

If when using the medical device, a power outage can create an unacceptable risk, the user and the installer must ensure that the medical device is connected to an appropriate power source such as an uninterruptable power supply.

3.5 Fixing the medical device to a non-removable support

Your medical device is not designed to be moved. To avoid accidentally dropping the device, we recommend that you find a suitable place to install it in your practice and attach it with the kit [F57811] supplied with it, to ensure that it cannot be dismantled or moved without using a tool.

4 Medical device description

4.1 Control unit

The control unit is compatible with:

1. The I-surge LED rotary technology via an I-surge LED surgical micromotor.
2. Newtron[®] ultrasonic piezoelectric technology via the surgical Piezotome[®] M+ handpiece or via the Newtron[®] dental handpiece and cord.

The piezo-ultrasonic technology common to both piezoelectric handpieces, called Newtron[®], is a patented system that automatically controls the frequency and power of tips in real time: the Cruise Control[®] system. Ultrasonic tips therefore maintain an even efficiency on the work area, regardless of the environment. Ultrasonic tips are active on hard tissue only, with minimal risks for soft tissues.

The micromotor can be used in combination with either the Piezotome[®] M+ handpiece, or with the Newtron[®] handpiece.

The touch-sensitive interface calls up all programmes stored in the memory and adjusts them before and after use.

4.2 Handpieces and I-Surge LED micromotor

Only SATELEC[®] handpieces and I-Surge LED micromotor can be connected to the medical device. For further information, refer to the documentation listed in chapter *Associated documentation* page 7.

4.3 Cords

The micromotor cord is only compatible with the I-Surge LED micromotor.

The Newtron[®] handpiece cord is only compatible with the Newtron[®] handpiece. It ensures irrigation circulation and is used to connect the medical device to the Newtron[®] handpiece.

The Piezotome[®] M+ cord is only compatible with the Piezotome[®] M+ handpiece. It ensures circulation and is used to connect the medical device to the Piezotome[®] M+ handpiece.

4.4 Inlets

The inlets at the back of the medical device provide suitable ventilation of the control unit and must be kept clear.

4.5 Switch

The mains switch is used to switch the medical device on (position I) or off (position O).

4.6 Mains connector

The mains connector with its earthing pin is used to connect the device to the electrical network via a disconnectable mains cord.












4.7 Fuse recess

The mains cord connector holds the recess for two mains fuses designed to protect the medical device in the event of overvoltage or an internal fault.

4.8 Pump





The medical device is fitted with two peristaltic pumps. A single-use pump is connected to the irrigation line.







5 User interface

Pictograms shared by the three modes	Meaning
	Start/stop purge
	Displays the flow in ml/min
	Irrigation ON
	Irrigation OFF
	Adjusts the irrigation flow rate.
	Handpiece lighting on
	Handpiece lighting off
	Progressive footswitch
	Smooth footswitch
	Warning. The type of warning is specified by a pictogram at the centre of the triangle.
	Saves any changes

5.1 I-surge mode

Before using the I-Surge mode, read the setting recommendations specified by the contra-angle or handpiece manufacturer

Pictograms for the I-Surge mode	Meaning
	Current contra-angle - press the pictogram to select a different contra-angle or personalise a specific contra-angle
	Revolutions per minute- adjust this value using the + and - pictograms on the screen
	Newton per centimetre - press the pictogram to adjust the torque of the force applied at the outlet of the contra-angle
	In P1, P2 and P3, adjust the number of revolutions per minute directly. In P4, adjust the torque in real time.







Pictograms for the I-Surge mode	Meaning
	Select programme P1 - Marking of the implant site
	Select programme P2 - Pilot drilling
	Select programme P3 - Boring/drilling
	Select programme P4 - Screwing
	Clockwise rotation
	Anti-clockwise rotation

I-Surge mode programme settings

Programme	Contra-angle ratio	Speed at end of tool	Torque at end of tool	Irrigation	Function
P1	20:1	1200 rpm	80 N.cm	80 ml/min	Implant site marking
P2	20:1	800 rpm	80 N.cm	100 ml/min	Pilot drilling
P3	20:1	15 rpm	20 N.cm	100 ml/min	Boring/drilling
P4	20:1	30 rpm	20 N.cm	0	Screwing

5.2 Piezotome® M+ mode

Before using an M+ tip, read the M+ tip kit user manuals and the Quick Reference guide for M+ tips [J87621].

Pictograms for the Piezotome mode	Meaning	Power
	Select programme D1 - Osteotomy, osteoplasty	Very powerful
	Select programme D2 - Osteotomy, osteoplasty	Powerful
	Select programme D3 - Osteotomy, osteoplasty	Average power
	Select programme D4 - Separation of soft tissue	Low power
	Used power per programme - press the + and - pictograms to adjust the value	
	Adjust the power	

Programme	Power	Irrigation
D1	3	60 ml/min

Programme	Power	Irrigation
D2	3	60 ml/min
D3	3	60 ml/min
D4	3	60 ml/min

Fine setting adjusts the power level for each programme between 1 and 5. The power value is a percentage of the maximum power at D1-5.







Programme	D1				
Fine setting level	1	2	3	4	5
Power value	82%	86%	90%	95%	100%
Frequency modulation	60 Hz				

Programme	D2				
Fine setting level	1	2	3	4	5
Power value	64%	68%	72%	74%	78%
Frequency modulation	60 Hz				

Programme	D3				
Fine setting level	1	2	3	4	5
Power value	44%	48%	52%	56%	60%
Frequency modulation	60 Hz				

Programme	D4				
Fine setting level	1	2	3	4	5
Power value	41%	46%	51%	56%	60%
Frequency modulation	30 Hz				


5.3 Newtron® mode

Pictograms for the Newtron mode	Meaning	Power
	Soft Mode	Low
	Medium Mode	Medium
	High Mode	High
	Boost Mode	Very high
	Used power - use the + and - pictograms to adjust the value	
	Adjust the power	

6 Medical device use

6.1 Preparing the medical device

6.1.1 Connecting the medical device to the electrical network

- Set the ON/OFF switch to 0 position (off).
- Connect the power cord to the device's mains connector.
- Connect the mains power cord to the mains socket equipped with an earthing pin.
- If necessary, connect the potential equalisation cable to the medical device's potential equalisation terminal marked with .

6.1.2 Footswitch

The practitioner can start up the medical device using the footswitch. Press the footswitch to automatically start the ultrasounds on Piezotome® M+ and Newtron® handpieces. The footswitch also controls the I-Surge LED micromotor.

The footswitch equipped with its cord cannot be disconnected. Its weight and antislip pad ensure good stability. The light function remains active for approx. 9 seconds after the footswitch is released.

6.1.3 Installing the footswitch

The footswitch must be positioned near the operator's feet and must always be within reach.

6.1.4 Connecting the I-Surge LED micromotor

Check for any traces of humidity on the I-Surge LED cord connector. Wipe them off if there are any.

Connect the cord connector to the socket, by aligning the red indexing points and by avoiding rotation movement. Place the I-Surge LED micromotor on the support.

6.1.5 Connecting the Piezotome® M+ handpiece

Check for any traces of humidity on the handpiece cord connector. Wipe them off if there are any.

Connect the cord connector to the socket, by aligning the red indexing points and by avoiding rotation movement. Place the handpiece on the support.

6.1.6 Connecting the Newtron® handpiece

Check for any traces of humidity on the handpiece cord connector. Wipe them off if there are any.

Connect the cord connector to the socket, by aligning the red indexing points and by avoiding rotation movement. Place the handpiece on the support.

6.1.7 Installing the irrigation system

Equipment: one sterile M+ irrigation line, 10 single-use sterile clips, one sterile perforator with cap, one bracket, one solution bag.


1. Insert the bracket into its holder
2. Open an irrigation line with its cassette
3. Slide the cassette into its support on the side of the medical device
4. Install the clips on the handpiece cord
5. Attach the irrigation line
6. Connect the end of the irrigation line to the water inlet on the handpiece
7. Install a bag on the bracket
8. Remove the cap from the perforator and insert the perforator into the solution bag.

| The medical device must always be used with SATELEC® M+ irrigation lines.

The medical device is not designed to administer medication and must only be used with bottles or bags containing no more than 1 litre of physiological saline solution or sterile water.

6.1.8 Purging the medical device

The medical device must be purged before and after use.




1. Immerse the irrigation line perforator in a container with distilled water.
2. Press  to rinse the irrigation line and the I-Surge LED micromotor or handpiece.
3. Minimum purge time is 1 minute.

6.2 Switching on the medical device

Place the I/O switch to I. The ImplantCenter™ M+ indicates that the device is switched on by an audible signal. The ImplantCenter™ M+ detects which handpiece cord is connected (Piezotome® or Newtron®) and the relevant menu will be displayed on the screen.

6.3 Selecting the operating mode


The ImplantCenter™ M+ can be used in any of the following modes:


-  I-Surge: with I-Surge LED micromotor and contra-angles or a straight handpiece
-  Piezotome: with Piezotome® M+ handpiece;
-  Newtron: with Newtron® handpiece.

Switch between modes by pressing the pictograms or using the footswitch as follows.

6.3.1 Common operations using the footswitch

The footswitch is used for common operations.






6.4 Connecting and disconnecting M+ accessories during use

- | Do not connect/disconnect the cords when the medical device is switched ON and the footswitch is pressed.
- | Do not tighten or loosen the contra-angles when the I-Surge LED micromotor is activated.
- | Do not tighten or loosen the tips when the handpiece is activated.

6.5 Switching off the medical device

| To avoid damaging the I-Surge LED micromotor, it must be rinsed with sterile water less than 30 minutes after use.

1. Remove the bag from the bracket.
2. Remove the irrigation line perforator from the bag.
3. Immerse the irrigation line perforator in a container with distilled water.
4. Press  and run the purge cycle for one minute to rinse the handpiece.
5. Remove the perforator from the container and purge the system until the medical device's irrigation line is completely empty.

6. Switch off the medical device (O/I switch to O).
7. Disconnect the I-Surge LED micromotor cord and the medical device's handpiece cord.
8. Remove the irrigation line clips and dispose of them in a special container for soiled medical equipment.
9. Disconnect the irrigation line from the I-Surge LED micromotor or handpiece and dispose of it in a special container for soiled medical equipment.
10. Disconnect the contra-angle or straight handpiece from the I-Surge LED micromotor and place it in the appropriate sterilisation box.
11. Rinse the micromotor and place it in the appropriate sterilisation box.
12. Detach the tip from the handpiece using the wrench and dispose of it in a special container for soiled medical equipment.
13. Remove the nozzle, optical guide and LED ring from the handpiece and place them in the appropriate sterilisation box.

This procedure is applicable for single-use sterile M+ tips. Non sterile M+ tips that can be used up to 5 times should be placed in the appropriate sterilisation box.






Refer to the cleaning, disinfection and sterilisation protocols for I-Surge LED micromotors, handpieces and tips for further information about the detailed sequence of pre-disinfection, cleaning, sterilisation operations for these types of medical device.

7 Configuring the medical device

In addition to normal use of the medical device, the following parameters can be configured:






- medical device volume
- screen brightness
- lighting timer
- factory settings

7.1 Adjusting the medical device volume

1. Switch on the medical device.
2. From the home screen, press  to open the menu.
3. Identify the pictogram that indicates volume .
4. Press  or  to adjust the volume.
5. Press  to save the settings and start the medical device.

Or press  to exit the settings screen without saving the changes.






7.2 Adjusting the screen brightness

1. Switch on the medical device;
2. From the home screen, press  to open the menu.
3. Identify the pictogram that indicates brightness .
4. Press  or  to adjust the screen brightness.
5. Press  to save the settings and start the medical device.

Or press  to exit the settings screen without saving the changes.




7.3 Adjusting the lighting timer

Press the footswitch to start the handpiece and release it to stop the handpiece. Once the handpiece is stopped, it will continue lighting up the treatment site during the time out period.


1. Switch on the medical device.
2. From the home screen, press  to open the menu.
3. Identify the pictogram that indicates lighting .
4. Press  or  to adjust the time out.
5. Press  to save the settings and start the medical device.

Or press  to exit the settings screen without saving the changes.

7.4 Recovering the medical device's factory settings

1. Switch on the medical device.
2. From the home screen, press  to open the menu.
3. Press  to recover all the factory parameters.
4. Press  to save the recovery and start the medical device.

Personalised contra-angles and changes made to programmes will be deleted from the medical device's memory.

Or press  to exit the settings screen without saving the changes.

7.4.1 I-Surge mode factory settings

Programme	Contra-angle ratio	Speed at end of tool	Torque at end of tool	Irrigation
P1	20:1	1200 rpm	80 N.cm	80 ml/min
P2	20:1	800 rpm	80 N.cm	100 ml/min
P3	20:1	15 rpm	20 N.cm	100 ml/min
P4	20:1	30 rpm	20 N.cm	0




7.4.2 Piezotome mode factory settings

	Programme	Power	Irrigation
Very powerful	D1	3	60 ml/min
Powerful	D2	3	60 ml/min
Average power	D3	3	60 ml/min
Low power	D4	3	60 ml/min

7.4.3 Newtron mode factory settings

	Programme	Power	Irrigation	Main functions
Green	Soft	P = 5	15 ml/min	Periodontics
Yellow	Medium	P = 5	15 ml/min	Endodontics
Blue	High	P = 5	15 ml/min	Scaling
Orange	Boost	P = 5	15 ml/min	Loosening

7.5 Personalising a contra-angle or handpiece

1. Switch on the medical device.
2. Press .
3. Press .
4. Enter the handpiece contra-angle ratio as specified by the manufacturer, using the keypad.
5. Rpm and Ncm values are recalculated automatically.
6. Press  to save the values.

8 Cleaning, disinfecting and sterilising

The instructions relating to cleaning, disinfection and sterilisation protocols for accessories provided by SATELEC® have been approved for each medical device and accessory. The applicable protocols are listed in the chapter *Associated documentation* page 7.

They can be downloaded at the following address:

www.satelec.com/documents

In all cases, the local regulations in force relating to the cleaning, disinfection and sterilisation protocols for accessories take precedence over the information provided by SATELEC®.

- | Do not use antiseptic products containing flammable substances.

8.1 Cleaning and disinfecting the medical device

The medical device must be OFF (switch on O) during cleaning and disinfecting procedures.

Avoid using cleaning and disinfection products that contain flammable agents. Otherwise, ensure that the product has completely evaporated or that there is no fuel left on the medical device and its accessories before switching it on.

- | Do not use an abrasive product to clean the medical device.

- | Never apply sprays directly to the medical device to clean it. Always spray the product onto a wipe, then clean the medical device.

Clean and disinfect the medical device's control unit and footswitch before and after use.

8.2 Cleaning, disinfecting and sterilising accessories

Refer to the cleaning, disinfection and sterilisation protocols for accessories listed in chapter *Associated documentation* page 7.

9 Monitoring and general maintenance of the medical device

Before and after use, check the medical device and its accessories entirely for any problems. This is necessary to detect any isolation fault or damage. If necessary, replace damaged parts.

Check that the air inlets on the control unit are clean to prevent any overheating.

| The clips holding the irrigation lines may cause wear of the handpiece cords. Check each handpiece cord individually before and after use.

10 Maintenance

Maintenance of the medical device essentially involves preventive maintenance operations, covering the following aspects:

- checking of accessories;
- everyday cleaning, disinfection and sterilisation ;
- cleaning.

10.1 Touch-sensitive screen messages

Use the touch-sensitive screen to configure the medical device. Depending on your actions, one or more of the following elements will be displayed.

10.1.1 Problem at start up

Symptom: the medical device beeps and displays the following pictogram:



Possible causes	Solutions
Internal error during start up: there is no communication between mother board and screen	<ul style="list-style-type: none">- Switch off the medical device (I/O switch to 0).- Wait 5 seconds before switching it back on.- Switch the device on (I/O switch to I) and don't touch any of the buttons while it is starting up.

10.1.2 Missing handpiece

Symptom: the medical device beeps when the user presses the footswitch and displays the following pictogram:



Possible causes	Solutions
Faulty connection between the handpiece cord and the medical device	<ul style="list-style-type: none">- Connect the handpiece cord to the medical device.- Switch off the medical device (I/O switch to 0).- Wait 5 seconds before switching it back on.- Switch the device on (I/O switch to I) and don't touch any of the buttons while it is starting up.

10.1.3 I-Surge LED micromotor missing

Symptom: the medical device beeps and displays the following pictogram:





Possible causes	Solutions
Faulty connection between the I-Surge LED micromotor cord and the medical device	<ul style="list-style-type: none"> - Connect the I-Surge LED micromotor cord to the medical device. - Switch off the medical device (I/O switch to O). - Wait 5 seconds before switching it back on. - Switch the device on (I/O switch to I) and don't touch any of the buttons while it is starting up.

10.1.4 Finding the software version

If you experience any problems with your medical device, the SATELEC® customer service team may ask you for the resident software version of your medical device.

To view the software version, proceed as follows:

1. Switch on the medical device.
2. From the home screen, press  to open the menu.
3. Note down the value displayed on the bottom of the screen.
4. Press  to start the medical device.

10.2 Identifying incorrect operation

In the event of incorrect operation, refer to the tables below to quickly identify and repair the non-complex parts of the medical device.

If the incorrect operation is not described in the tables below, please contact your supplier or the customer service team at SATELEC®.

Do not use the medical device if it appears to be damaged or faulty. Isolate the medical device and make sure that it cannot be used.

10.2.1 Not working

Symptoms: the screen is off, the medical device failed to start

Possible causes	Solutions
Faulty connection between the mains cord and the medical device	Connect the mains cord to the medical device
Faulty connection between the mains cord and the electrical wall socket	Connect the mains cord to the electrical wall socket
Switch in position O	Set the switch to position I
Mains fuses in the mains connector not working	Replace the mains fuses with fuses of the same type and rating
Internal fuse not working	Return the device to the SATELEC® customer service team
No electrical current	Contact your electrician
Transmission fault	Switch off the medical device then switch it on again Return the medical device to the SATELEC® customer service team

The medical device also has an internal fuse that cannot be accessed by the user.

10.2.2 The power is not as expected

Symptoms: the tip does not vibrate at the expected frequency, the treatment is not progressing as normal and is taking longer or at a standstill.

Possible causes	Solutions
Worn or bent tip	Replace the tip

Possible causes	Solutions
Power setting incorrect	Adjusts the power
Angle of approach incorrect or inadequate pressure on the clinical site	Refer to the tip user manual indicated chapter <i>Associated documentation</i> page 7
Moisture on the I-Surge LED micromotor or handpiece cord connector	Dry the electrical contacts

10.2.3 No spray or very little amount of spray

Symptoms: when the device is in use, irrigation is not working and no spray comes out of the contra-angle or the tip

Possible causes	Solutions
Irrigation solution bag empty	Replace the irrigation bag with a new one
Irrigation deactivated	Start irrigation using the footswitch or touch-sensitive screen
Irrigation line blocked	Replace the irrigation line
Irrigation line pinched	Check the length of the irrigation line
Incorrect irrigation adjustment	Adjust the irrigation flow rate

10.2.4 Light not working or too dim

Symptoms: the handpiece is not lighting up the clinical site

Possible causes	Solutions
LED ring missing in the handpiece nozzle	Install the LED ring
Faulty LED ring	Replace the LED ring
Cracks or fine cracks on the LED ring	Replace the LED ring
LED ring contacts faulty	Replace the LED ring
LED ring poles reversed	Install the LED ring with the correct polarity
Optical guide damaged	Replace the optical guide
Handpiece cord connector faulty	Dry the cord and medical device connectors Replace the handpiece and cord
Other	Return the medical device to the SATELEC [®] customer service team

10.2.5 Water leakage

Symptoms: water is leaking from one of the following places:

- irrigation cassette
- handpiece

Possible causes	Solutions
The handpiece irrigation inlet is damaged	Replace the handpiece and cord Return the handpiece to the SATELEC [®] customer service team
Pipe ruptured in the irrigation line cassette	Replace the irrigation line
Irrigation cassette not working	Replace the irrigation line

10.2.6 Ultrasounds not working

Symptoms: the tip does not vibrate, vibration cannot be heard.

Possible causes	Solutions
Tip loose	Fasten the tip using the wrench
Connector with the medical device faulty	Clean the cord contacts
Handpiece cord wire(s) cut	Return the handpiece and cord to the SATELEC [®] customer service team to have them replaced

10.3 Corrective maintenance

In the event of faulty operation, the following corrective maintenance actions may be performed by the user.

10.3.1 Replacing the fuses

The medical device is protected by two fuses in the mains connector.

To replace the fuses, perform the following operations:

- stop the medical device (position O);
- disconnect the mains cord from the electrical network;
- disconnect the mains cord from the mains connector;
- insert the tip of a flathead screwdriver into the notch on top of the fuse holder to release it;
- remove the used fuses;
- replace the used fuses with fuses of the same type and same rating;
- place the fuse holder in its recess by pushing it until you hear a click that confirms it is in the correct position;
- connect the mains cord to the connector;
- connect the mains cord to the electrical network;

11 Technical specifications of the medical device

11.1 Identification

Manufacturer	SATELEC®
Name of the medical device	IMPLANTCENTER™ M+

11.2 Control unit

Width (in mm)	472.9
Height (in mm)	149.5 - 471.1 with bracket
Depth (in mm)	339.9
Weight	5 kg without accessories

Ingress protection rating: IPX0

11.3 Generator

Supply voltage	100 - 230 VAC
Power supply frequency	50 Hz / 60 Hz
Power consumption	250 VA to 230 VAC
No-load voltage	250 V
No-load voltage mode Newtron®	150 VA
Minimum output frequency	28 kHz
Ultrasonic operating mode: Piezotome® and Newtron®	Intermittent operation 10 minutes ON / 5 minutes OFF
I-Surge operating mode	Intermittent operation 20 s ON / 30 s OFF at 2 N.cm
I-Surge LED micromotor speed	100 to 40,000 rpm
I-Surge LED micromotor torque	Max. 6 N.cm
Electrical rating	1
Electrical safety class, leakage current types	LF
Internal fuse not accessible to the user	1 fuse - F1: 5 mm x 20 mm - 10 AT / 250 VAC
Fuse (mains connector)	2 fuses - 5 mm x 20 mm - 2 AT for 100 to 230 VAC

Thermal protection against very high temperatures of the I-Surge LED micromotor.

11.4 Screen

Width (in mm)	115
Height (in mm)	86

11.5 Length of cords

Newtron® handpiece cord (in mm)	2900
Piezotome® M+ handpiece cord (in mm)	2900
Footswitch cord (in mm)	2900

11.6 I-Surge LED micromotor

Length (in mm)	90,6
Maximum diameter (in mm)	23.2
Weight (in g)	119 without cable
Coupling	as per ISO standard 3964

11.7 Irrigation

Nominal water output flow at the handpiece tip of the ImplantCenter™ M+ in ml/min	10 to 120 ml
Nominal water output flow at purge in ml/min	120 ml per min
Piezotome® mode settings	10 ml/min pitch
Newtron® mode settings	1 ml/min pitch

11.8 Footswitch

Width (in mm)	311
Height (in mm)	181
Depth (in mm)	209
Weight	3.5 kg

Ingress protection rating: IPX8

11.9 Environmental characteristics

Operating temperature	10 to 30° C
Storage temperature	0 to 50° C
Operating humidity	30 to 75%
Maximum storage humidity	100%, including condensation
Atmospheric pressure	500 hPa to 1060 hPa
Altitude	2000 metres maximum

11.10 Environmental restrictions

Usage premises	The medical device is designed to be used in an operating theatre or in premises suitable for the procedure involved.
Use in gas-filled atmosphere	The medical device is not designed to be used in an AP or APG gas-filled atmosphere.
Immersion	The medical device must not be immersed.

11.11 Main performance characteristics

- Ultrasonic vibrations of the surgical tip fitted to the end of the surgical handpiece.
- Vibration frequency 28 kHz or greater.
- Modulation frequency

12 Regulations and standards










12.1 Official texts











This medical device complies with the essential requirements of European Directive 93/42/EEC. This equipment is designed and developed in compliance with the Electrical Safety standard IEC60601-1 in force. It was designed and manufactured in accordance with an EN ISO 13485-certified quality assurance system. This documentation complies with European regulation No.207/2012.

12.2 Medical class of the device

Class of medical device: IIb according to 92/42/EEC directive

12.3 Symbols

Symbols	Meaning
	Authorised humidity range
	Authorised temperature range
	Authorised pressure range
	Quantity (1)
	Keep away from humidity
	Fragile
	Not to be used on patients with implantable medical devices
	Refer to the accompanying documentation
	Consult the User Manual

Symbols	Meaning
 Electronic user informations	Accompanying documentation in electronic format
	LF type
I	Class 1
	Alternating voltage
	Sterilisation at 134 °C in an autoclave
	Sterilisation at 132 °C in an autoclave
	Washer disinfector for thermal disinfection
	EC marking
	Do not dispose of as household waste
YYYY 	Year of manufacture
	Footswitch
O	Device OFF
I	Device ON
IPX8	IP: ingress protection ratings procured by a range X: no ingress of protection rating claim against the penetration of solids 8: protected against the effects of continuous immersion in water

12.4 Manufacturer identification

SATELEC
 A Company of ACTEON Group
 17, avenue Gustave Eiffel
 BP 30216
 33708 MERIGNAC cedex
 FRANCE
 Tel. +33 (0) 556.34.06.07
 Fax. +33 (0) 556.34.92.92
 E.mail: satelec@acteongroup.com.
 www.acteongroup.com



12.5 Branch addresses

U.S.A. & Canada

ACTEON North America
124 Gaither Drive, Suite 140
Mount Laurel, NJ 08054 - USA
Tel. +1 856 222 9988
Fax. +1 856 222 4726
info@us.acteongroup.com

GERMANY

ACTEON GERMANY GmbH
Industriestrasse 9 - 40822 METTMANN - GERMANY
Tel. +49 21 04 95 65 10
Fax. +49 21 04 95 65 11
info@de.acteongroup.com

SPAIN

ACTEON MEDICO-DENTAL IBERICA, S.A.U.
Avda Principal n° 11 H
Poligono Industrial Can Clapers
08181 SENTMENAT (BARCELONA) - SPAIN
Tel. +34 93 715 45 20
Fax. +34 93 715 32 29
info@es.acteongroup.com

U.K.

ACTEON UK
Unit 1B - Steel Close - Eaton Socon, St Neots
CAMBS PE19 8TT - UK
Tel. +44 1480 477 307
Fax. +44 1480 477 381
info@uk.acteongroup.com

MIDDLE EAST

ACTEON MIDDLE EAST
Numan Center - 2nd Floor N° 205 - Gardens Street
PO Box 5746 - 11953 AMMAN - JORDAN
Tel. +962 6 553 4401
Fax. +962 6 553 7833
info@me.acteongroup.com

CHINA

ACTEON CHINA
Office 401 - 12 Xinyuanxili Zhong Street -
Chaoyang District - BEIJING 100027 - CHINA
Tel. +86 10 646 570 11/2/3
Fax. +86 10 646 580 15
beijing@cn.acteongroup.com

THAILAND

ACTEON (THAILAND) LTD
23/45 Sorachai Building 16th floor - Sukumvit 63

Road, Klongton Nua - Wattana, BANGKOK 10110
- THAILAND

Tel. +66 2 714 3295
Fax. +66 2 714 3296
info@th.acteongroup.com

INDIA

ACTEON INDIA
B-94, GIDC Electronic Estate - Sector 25 -
GANDHINAGAR 382028 Gujarat - INDIA
Tel. +91 79 2328 7473
Fax. +91 79 2328 7480
info@in.acteongroup.com

LATIN AMERICA

ACTEON LATINA AMERICA
Bogotá - COLOMBIA
Mobile: +57 312 377 8209
carlosandres.vera@es.acteongroup.com

RUSSIA

ACTEON RUSSIA
Valdajski Proezd 16 - office 243
125445 Moscow - RUSSIA
Tel./Fax. +7 499 76 71 316
sergey.koblov@ru.acteongroup.com

AUSTRALIA/NEW ZEALAND

ACTEON AUSTRALIA/NEW ZEALAND
Suite 119, 30-40 Harcourt Parade
Rosebery NSW 2018
Australia
Tel. +612 9669 2292
Fax. +612 9669 2204
sandy.junior@au.acteongroup.com

TAIWAN

ACTEON TAIWAN
11F., No.1, Songzhi Rd.
Xinyi Dist., Taipei City 11047
TAIWAN (R.O.C.)
+ 886 2 8729 2103
tina.chu@tw.acteongroup.com

13 Disposal and recycling

As an item of Electrical and Electronic Equipment, the medical device must be disposed of via a specialist collection, removal, recycling or destruction channel. This applies in particular to the European market, in reference to Directive 2002/96/EC dated 27/01/2003.

When your medical device has reached the end of its service life, contact your nearest equipment dealer, or ACTEON GROUP head office or one of the company branches to find out how to proceed. The relevant contact details are given in chapter *Branch addresses* page 35.



14 Index

A

approved dealers	10
attachment	11
attachment system	
risk of fall	12

D

Date	9
------------	---

E

electrical safety	33
earthing pin	13
environment	32
gas-filled atmosphere	32
pressure	32
temperature	32

F

first inclusion of EC marking	9
-------------------------------------	---

H

humidity	17
----------------	----

I

I-Surge LED	7
incorrect operation	28
damage	25
fault	25

fuses	13, 28, 30
water leakage	29
indexing points	17

M

Manufacturer	31
medical class	33
medical device	
air inlets	13
Cruise Control®	13
footswitch	17
light function	17
mains connector	13, 30
peristaltic pump	13
settings	22
switch	13, 28

R

recycling	37
disposal	37
repair	10

S

surgery	9
---------------	---

T

tip	29
sterile	10

U

user documentation	
download	7
electronic user instructions	7
European regulations	33

instructions for the entire range of M+	
ultrasonic generators	5
paper documentation	7
Quick Clean	7
Quick Start	7
update	9
user manual	7

15 Glossary

with the device, or instructions for use available through a website [European regulation No.207/2012]

A

accessory

includes micromotors, cords, hand-pieces, nozzles, LED rings, optical guides, tips, sterilisation boxes, and silicone supports used in conjunction with the control units

B

bag

container filled with physiological saline solution or sterile water. Can designate either a flexible container or rigid container like a bottle.

F

fixed installed medical devices

devices and their accessories which are intended to be installed, fastened or otherwise secured at a special location in a healthcare facility so that they cannot be moved from this location or detached without using tools or apparatus, and which are not specifically intended to be used within a mobile healthcare facility [European regulation No.207/2012]

I

Instructions for use

information provided by the manufacturer to inform the user of the device of its safe and proper use, of its intended performances and of any precautions to be taken [European regulation No.207/2012]

instructions for use in electronic form

instructions for use displayed in electronic form by the device, contained in portable electronic storage media supplied by the manufacturer together

P

professional users

persons using the medical device in the course of their work and in the framework of a professional healthcare activity [European regulation No.207/2012]



Ref - J27251 • V2 • 01/2016 • NO25EN010B

