Motorsystem MD 10

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Nr. 31820

www.nouvag.com

NOUVAG®
Dental and Medical Equipment

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Congratulations on your purchase of a NOUVAG AG product. Thank you for the confidence shown in our products.

Please consult the instruction manual for the use and maintenance of the device in order to ensure that it will function properly and efficiently for many years.

You will find the conformity statement and list of authorized service representatives attached.

- **Before operating, please read instructions carefully!**

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TABLE OF CONTENTS

1	Description	2
1.1	Operations and functions	2
1.2	Technical Data	2
1.3	Operating Environment	2
1.4	Transport and Storage Environment	2
1.5	Symbols	2
2	Safety measures	3
3	Supplied equipment	3
4	System Assembly	4
5	Operation	4
5.1	Set-up	4
5.2	Tubing assembly-Internal irrigation	5
5.3	Tubing assembly-External irrigation	5
5.4	Tubing assembly-combined	5
6	Operation	6
6.1	Starting the MD 10	6
6.2	Keyboard	6
6.3	Setting the revolutions-	7
6.4	Setting the output of the irrigation pump	7
6.5	Operation with vario-pedal	7
6.6	Setting the torque limiter AL (Automatic Limiter)	8
6.7	Setting the torque limiter AS (Automatic Stopper)	8
6.8	Memory	8
7	Disinfection, Cleaning and Sterilization	9
7.1	Control unit and foot control	9
7.2	Electronicmotor	9
7.3	Tubing set Nr. 1706	11
7.4	Y-Connector	11
8	Troubleshooting	11
9	Replacement parts and numbers	12
10	Technical Safety Inspections	12
11	Disposal	12

1 Description

1.1 Operations and functions

The MD10 is used in Implantology and microsurgery. The unit is designed for drilling, screwing and sawing. A built-in irrigation system minimizes the heat generation of the rotating instruments in order to reduce infections that can cause tissue damage.

1.2 Technical Data

Voltage: changeable: 230 V~/ 50-60 Hz or 115 V~/ 50-60 Hz
 Power: 130 VA
 Application of part: Type BF
 Protection level: Class 1
 Motor speed: 500 – 40,000 rpm
 Max. torque: 10 Ncm
 Dimensions (WxHxD): 120 x 190 x 250 mm
 Net weight: 3 kg

1.3 Operating Environment

Relative Humidity: max. 80%
 Temperature: 10 to 40°C
 Pressure: 800 to 1060 hPa

1.4 Transport and Storage Environment

Relative Humidity: max. 90%
 Temperature: 0 to 60°C
 Pressure: 700 to 1060 hPa

1.5 Symbols

: Attention accompanying documents

: suitable for thermal disinfection

: Application of part type BF

IP68

: Watertight foot control

CE 1275 : Conforms to EU standards

KB 1min on / 3min off : Electronic motor short operation:
1min on 3min off

C_{US} : Certified by Canadian Standards Association (CSA) for Canada and USA

: Old electrical and electronic equipment must be disposed separately and may not be included in regular domestic waste.

: Autoclavable at max. 135°C

2 Safety measures

Your safety, the safety of your team, and, it goes without saying, the safety of your patients is for Nouvag AG the first priority. It is therefore vital that the following measures be strictly observed:

Important:

- The MD 10 should be operated by qualified personnel only!
- Responsibility for the use of sub-units, accessories, parts or assemblies from other manufacturers rests solely with the user!
- Repairs are to be carried out by authorized NOUVAG service technicians only!
- Nouvag AG cannot be held liable for any malfunction of the MD 10, or performance failure and/or its designed or desired utility, nor can Nouvag AG be held liable for any injuries to persons or animals, in any case when the MD 10 is miss-used or not operated, applied or maintained in strict accordance with the user/owner instructions set out in the operating manual. In the event of any doubt or question, the user is to contact Nouvag AG or its lawful representative for clarification or assistance!
- For care please use the Nou-Clean Spray. Use of not equivalent care products could cause to malfunction and loss of warranty.

In Operation:

- The unit will be not supplied in a sterile condition!
- Do not attach the handpiece and contra angle while the motor is running!
- To avoid danger, do not grab running burr or drill!
- To avoid damage of the instruments, do not handle the clamping mechanism, while the handpiece or contra angle are in operation!

3 Supplied equipment


Art. No.	Description	Quantity
2150	Control unit MD 10	1 unit
1866	Vario footpedal; IP68; electr.....	1 unit
2063	Electronic motor 31 ESS	1 unit
6024	Hose set, sterile, 3 m, one-way.....	1 unit
1696	Coolant; 0.9% NaCl, 1 l.....	1 unit
1770	Support for coolant.....	1 unit
1942	Add-on to Nou-Clean Spray; for electronic motor 31 ESS	1 unit
1958	Add-on to Nou-Clean Spray; for instruments.....	1 unit

Following articles cannot be included within the delivery of the Control unit, due to law restrictions in delivery of dangerous goods:

1984	Nou-Clean; care and maintenance spray	1 unit
------------	---	--------

4 System Assembly

- | | |
|---------------------------------|----------------------------------|
| 1. Control unit | 12. Pump |
| 2. Keyboard panel | 13. Mains cable |
| 3. Mains switch "ON / OFF" | 14. Foot control (IP68) |
| 4. Pedal socket "FOOT SWITCH" | 15. Switch for Pump and Forw/Rev |
| 5. Motor socket "MOTOR" | 16. Foot plate |
| 6. Electronicmotor | 17. Tube clamp |
| 7. Handpiece support | 18. Drip chamber |
| 8. Contra angle (not included) | 19. Air valve |
| 9. Tubing set | 20. Mains cable socket |
| 10. Cooling fluid flask | 21. Voltage switch |
| 11. Hanger rod for liquid flask | |


5 Operation

5.1 Set-up

- 1) Insert solution hanger in support on panel of control unit.
- 2) Insert electronicmotor plug into motor socket "MOTOR".
- 3) Insert pedal plug into pedal socket, "FOOT SWITCH".
- 4) Attach contra angle to electronicmotor.

→ see row of pictures right:

- 5) Loosen and lift screw latch from pump lid.
- 6) Open lid forward.
- 7) Remove demonstration tube.
- 8) Insert tube set.
- 9) Close pump lid, lower screw latch and close.


Insure that „pressure“ tube runs to irrigation needle, and „suction“ tube runs to solution bottle!


- 10) Insert tube needle and hang cooling flask on hanger rod.
- 11) Open tube clamp.
- 12) Open air valve of the drip chamber.
- 13) Connect the mains cable into the electrical outlet.

Check for correct operation voltage and verify mains voltage output!


5.2 Tubing assembly-Internal irrigation

- 1) Insert internal irrigation needle (1) into the end of the tubing.
- 2) Attach clip (2 & 3) to irrigation tubing.
- 3) Insert irrigation needle into contra angle.
- 4) Attach clip to contra angle and to motor cable.
- 5) When necessary, attach additional clips to motor cable.


5.3 Tubing assembly-External irrigation

- 1) Attach clip (1 & 2) to irrigation tubing.
- 2) Connect piece of tubing to irrigation needle (3).
- 3) Attach clip to contra angle and to motor cable.
- 4) When necessary, attach additional clips to motor cable.

5.4 Tubing assembly-combined

- 1) Attach Y-connector (1) to the end of the tubing.
- 2) Attach two 16cm pieces of tubing (2) to Y-connector.
- 3) For internal irrigation, connect first tubing piece to irrigation needle (3).
- 4) Attach clips (4 & 5) to irrigation tubing.
- 5) Attach irrigation tubing with clip to contra angle and to motor cable.
- 6) Insert irrigation needle into contra angle.
- 7) Connect second tubing piece with external irrigation.
- 8) When necessary, attach additional clips to motor cable.


6 Operation

6.1 Starting the MD 10

Press green "ON/OFF" switch.

6.2 Keyboard

- A) **Display:** Indication of the revolutions (speed), torque or pump speed.
- B) **Multiplying key "1:5"** (for contra angle 1:5):
Speed: 2,500 – 200,000 rpm Torque: 1, 2 Ncm
- C) **Multiplying key "1:1"** (for contra angle 1:1):
Speed: 500 – 40,000 rpm Torque: 1, 2, 3, 4, 5, 6, 7, 8, 9, 10 Ncm
- D) **Reduction "16:1"** (for contra angle 16:1):
Speed: 30 – 2,500 rpm Torque: 5, 7, 10, 12, 15, 17, 20, 25, 30, 35 Ncm
- E) **Reduction "32:1"** (for contra angle 32:1):
Speed: 15 – 1,250 rpm Torque: 10, 15, 20, 25, 30, 35, 40, 45, 50, 55 Ncm
- F) **Fixed key "SAW"** (for Nouvag saws):
Speed: 14,000 rpm Torque: 10 Ncm
- G) **Fixed key "MUCO"** (for Nouvag Mucotome):
Speed: 8,000 rpm Torque: 10 Ncm
- H) **Key "ATC SET":**
Press key "ATC SET" to adjust the torque limiter AL or AS.
- I) **Key "PUMP SET":**
Press key "PUMP SET" to adjust the pump speed.

The green light on key B) to I) illuminates, when the key is activated.

- J) **Taste "DOWN":**
To decrease the speed, torque and pump.
- K) **Taste "UP":**
To increase the speed, torque and pump.
- L) **Light "MOTOR ON":**
The green light indicates that the electronicmotor is switched on.
- M) **Light "REVERSE MODE":**
The green light indicates that the electronicmotor operates in reverse and a warning tone will sound.

6.3 Setting the revolutions-

- 1) Press the appropriate reduction or multiplying key (for example key "32:1" for contra angle 32:1).
- 2) Press keys "UP" or "DOWN" to select the desired speed.

With the keys "SAW" (for Nouvag-saws) and "MUCO" (for Nouvag-Mucotome) the corresponding fixed speed and torque values can be adjusted.

The reduction or multiplication on the keyboard panel has to be according to the reduction or multiplication of the contra angle, otherwise the device will indicate an incorrect revolution.

6.4 Setting the output of the irrigation pump

The output of the irrigation pump can be selected from 1 to 10 levels.

- 1) Press key "PUMP SET" and in the display appears for example, pump selection [P 6].
- 2) Press keys "UP" and "DOWN" to select the desired pump speed.
- 3) Press key "PUMP SET" or wait 3 seconds until the revolutions appear automatically in the display.

6.5 Operation with vario-pedal

Key "PUMP":

To turn the pump on and off (see light "PUMP SET")

Key "FORW / REV":

To adjust the direction of rotation of the motor (see light "REVERSE MODE")

Foot plate:

Push foot plate down to start pump action and to vary the motor speed.

Foot plate ...	Motor:	Pump:
... not pressed	Motor off	Pump off
... lightly pressed	Motor runs slowly	Pump on, if "PUMP" is pressed (speed as indicated by setting on control unit)
... fully pressed	Motor speed max.(as indicated by setting on control unit)	Pump on, if "PUMP" is pressed (speed as indicated by setting on control unit)

Caution: for safety use the unit with foot control only!

6.6 Setting the torque limiter AL (Automatic Limiter)

The automatic limiter (AL) limits the torque level applied to the instrument. The AL-modus is used for drilling into the bone, in order to allow the electronic motor to develop enough power, so the instrument can have a good drive in hard bones, as well. The speed in the instrument is held constantly until the selected setting is reached. If the load in the instrument over the adjusted limiter increases, the speed will drop to zero. The power in the instrument remains. If the load in the instrument over the adjusted limiter decreases, the speed will increase again.

All torque values are shown in Ncm (for example a setting of [AL 20] → indicates that the instrument can be loaded to maximum of max. 20 Ncm before rotations stops).

Adjusting the AL:

- 1) Press the "ATC" key. In the display appears the torque mode set (for example [AL 20]).
- 2) Pressing the keys "UP" and "DOWN" will increase and decrease the torque level respectively.
- 3) Press "ATC SET" or wait for 3 seconds until the speed appears automatically in the display.

If the display says [AS.....] :

Press key "DOWN" until [AL.....] appears in the display.

When the AL-Modus is activated, the green light on the "ATC SET" key does not illuminate.

6.7 Setting the torque limiter AS (Automatic Stopper)

The function AS limits the torque in the instrument. The AS operates like a torque wrench for tightening implants and abutments. The electronic motor stops immediately when the preselected torque is achieved and in the display appears the selected torque, for example [AS 40]. The electronic motor does not generate power anymore. To restart the electronic motor, release briefly the foot plate and press it again.

All torque values are shown in Ncm (for example a setting of [AS 20] → indicates that the instrument can be loaded to maximum of max. 20 Ncm.

The function AS is activated in the reductions 16:1 and 32:1 only, and up to 150 rpm and 75 rpm respectively! For other reduction or multiplication use function AL.

Setting the AS:

- 1) Press reduction key "16:1" or "32:1".
- 2) Adjust the speed between 30 and 150 rpm and 15 and 75 rpm respectively.
- 3) Press key "ATC SET" . In the display appears the torque selection (for example [AS 20]).
- 4) Press keys "UP" and "DOWN" to select the desired torque values.
- 5) Press the key "ATC SET" or wait for 3 seconds until the revolutions appear automatically in the display.

If the display says [AL.....]:

Press key "UP" until [AS.....] appears in the display.

When the AS-Modus is activated, the green light on the "ATC SET" illuminates.

6.8 Memory

Once programmed, the preselected speed, torque and value setting will remain stored in memory, even when the motor is turned off.

7 Disinfection, Cleaning and Sterilization

Please pay attention to the following important points for the maintenance of the material:


- **Clean, disinfect and sterilize the device and components after each use.**
- **Do not use dissolving agents for cleaning.**
- **Autoclave material in transparent packaging.**
- **Do not fill the sterilization bag more than 80%.**
- **Autoclave material at maximum 135°C.**
- **Sterilized material should be stored and tagged with sterilization date,**

7.1 Control unit and foot control

The control unit and the foot control do not come in contact with patients. Clean the surface of the unit only, do not use harsh cleaners or solvents for cleaning. Use 80% ethyl alcohol or microbiologically effective disinfectants. The control unit face is sealed and washable.

7.2 Electronicmotor

Restrictions of treatment	Frequent reprocessing hardly affects the Electronicmotor. The end of the product lifetime is usually determined by damages and wear and tear caused by usage.
----------------------------------	---

INSTRUCTIONS	
Usage location:	Wipe dirty surfaces with an expendable cloth/ tissue paper.
Storage and transport:	No special requirements. Avoid lengthy waiting periods before preparation because of danger of drying and corrosion.
Cleaning preparations:	Remove any surface soiling of the Electronicmotor with a disposable/paper towel. Unscrew motor cover and remove it with the cable. 
Automatic cleaning and disinfecting:	Equipment: Cleaning/disinfection equipment with special load carrier, which ensures the connection of the motors to the cleaning/disinfecting equipment and the irrigation of the channels. Wash the Electronicmotor from the front side. Use neutral cleaners as cleansing agents. 1) Load the Electronicmotor in the load carrier (irrigation of the channels must be ensured). 2) Load Motor cover with cable in basket. 3) Set cleaning cycle with effectual rinsing and cleaning. Carry out rinse cycle with demineralised water. 4) Rinse with demineralised water at 93°C for 10 minutes for thermal disinfecting. 5) Check whether there is still visible dirt in the gaps when taking the Electronicmotor and Motor cover with cable out. If necessary repeat cycle or clean manually.

Manual cleaning:	Equipment: Neutral cleaning agent, soft brush, running water. Procedure: 1) Rinse and brush surface soiling from the Electronicmotor and Motor cover with cable. 2) Apply cleaning agent with brush on all surfaces as well as gaps. 3) Rinse the Electronicmotor and motor cover with cable thoroughly in running water.
Manual disinfecting:	Spray the Electronicmotor and motor cover with disinfection spray on ethyl alcohol basis.
Drying:	Dry the Electronicmotor and motor cover with cable manually or in a drying chamber, if the drying/disinfecting equipment does not provide any drying program.
Inspection and maintenance:	Check for visual damage, corrosion and wear and spray the Electronicmotor with NOU-CLEAN-SPRAY thoroughly and briefly wipe with a moistened cloth (see illustration and instructions on the spray). Screw motor cover with cable on the Electronicmotor after spraying. 
Packaging:	Single: Pack the Electronicmotor in an individual sterile supply package. The bag must be large enough so that the seal is not under tension. Sets: Arrange the Electronicmotors into the foreseen trays or put onto an all-purpose sterilisation tray.
Sterilisation:	Autoclave in vacuum autoclaves at max. 135°C for at least 5 minutes*. Do not exceed the maximum load of the steriliser, if several instruments are sterilised during one sterilisation cycle. Observe a drying phase for autoclaves without post-vacuum. Allow the Electronicmotor to dry in the bag with the paper side upwards. *The temperature hold times must comply with the local guidelines and standards.
Storage:	No special requirements. The sterilisation indicator and sterilisation date should be marked on the packing, if the sterilised Electronicmotor is not used immediately after sterilisation.

The above-mentioned instructions were validated as suitable from NOUVAG AG. The person accomplishing the treatment is responsible that the treatment that actually took place leads to the desired results in the treatment facility with the available equipment, materials, and personnel. To ensure this it is normally necessary to validate and supervise the process. Each deviation from the instructions provided should be evaluated carefully with regard to effects and negative consequences.

7.3 Tubing set Nr. 1706

Disposable tubing set Nr. 1706 can not be sterilized!

7.4 Y-Connector

After use, flush the saline residue out of the Y-Connector with distilled water for approximately 20 seconds with distilled water!

- 1) Flush the Y-Connector with connected contra angle with distilled water for approximately 20 seconds.
- 2) Remove the tubings from the Y-Connector.
- 3) Pack the Y-Connector in a transparent bag and sterilize the pack in an autoclave at **max. 135°C**.

8 Troubleshooting

Problem:	Cause:	Solution:	Refer: operating Manual:
Device does not operate at all	Control unit not turned on	Turn main switch "ON/OFF" ON	4 System Assembly
	Wrong operating voltage	Check mains voltage output	1.2 Technical Data
	Mains not connected	Connect mains to control unit	4 System Assembly
No coolant in the instrument	Irrigation pump not turned ON	Turn the irrigation pump ON	6.5 Operation with vario-pedal
	Tubing set not mounted properly	Mount tubing properly (pay attention to direction!)	5.1 Set-up
	Tubing set is glued together/coated with deposit	Exchange tubing	5 Operation
	Cooling fluid flask not ventilated	Open air filter on drip chamber	4 System Assembly
	Tubing set is dripping	Exchange tubing set	5.1 Set-up
[HOT] appears in display	Electronic motor is getting too hot	Turn off the unit, wait for 10 minutes, turn on the unit again	
Foot control fails to function	Foot control not connected	Plug foot control into control unit	5.1 Set-up

In the event a problem cannot be solved, contact dealer or authorized service technicians, as listed on the last page of your Operating Instruction Manual.

9 Replacement parts and numbers

Accessories	Art.-Nr
Clip Set motor cable (10 pcs.)	1873
Clip Set contra angle (3 pcs.)	1881
Tubing (16cm) for Y-Connector	1773
Disposable tubing set, standard, sterile	1706
Y-Connetor	1777
Cooling fluid flask, 0.9% NaCl, 1l	1707
Nou-Clean; Cleaning spray	1984
Spray nozzle attachment E-coupling, for surgery instruments.....	1958
Spray nozzle attachment, for electronicmotor 31 ESS	1942

For ordering other parts, our customer service personnel are pleased to assist

10 Technical Safety Inspections

Various countries require technical safety inspections of medical devices in regulations. The technical safety inspection is a legal periodic safety check required for operators of medical devices. The objective of this measure is the timely identification of faults and risks for patients, users or third parties.

The Inspection interval for the MD 10 was stipulated on **2 years** (based on Appendix 1 of the MPBetreibV) [German Medical Devices Operator Ordinance]:

NOUVAG AG offers technical safety inspections for clients. Please contact our technical customer service.

11 Disposal

Disposal of device, components and accessories must strictly conform to local laws and regulations as set out by the relevant authorities.

With regard to the preservation of the environment old equipment may be returned to the distributor or manufacturer.

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Electromagnetic compatibility (EMC)

Remark:

The **Product** subsequently referred to herein always denotes the MD10.

Changes or modifications to this product not expressly approved by the manufacturer may result in increased emissions or decreased immunity performance of the product and could cause EMC issues with this or other equipment. This product is designed and tested to comply with applicable regulations regarding EMC and shall be installed and put into service according to the EMC information stated as follows.

WARNING

Use of portable phones or other radio frequency (RF) emitting equipment near the product may cause unexpected or adverse operation.

WARNING

The product shall not be used adjacent to, or stacked with, other equipment. If adjacent or stacked use is necessary, the product shall be tested to verify normal operation in the configuration in which it is being used.

Compliant Cables and Accessories

WARNING

The use of accessories, transducers and cables other than those specified may result in increased emissions or decreased immunity performance of the product.

The table below lists cables, transducers, and other applicable accessories for which the manufacturer claims EMC compliance.

NOTE: Any supplied accessories that do not affect EMC compliance are not listed.

Description	Length max.
Electronic motor 31ESS	2.9m
Vario Footpedal IP68	2.9m

Guidance and manufacturer's declaration – electromagnetic emissions

The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration – electromagnetic immunity

The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.

Immunity tests	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 8 kV air	+/- 6 kV contact +/- 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/- 1 kV differential mode +/- 2 kV common mode	+/- 1 kV differential mode +/- 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5 % U_T (> 95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (> 95 % dip in U_T) for 5 sec	< 5 % U_T (> 95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (> 95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Product requires continued operation during power mains interruptions, it is recommended that the Product be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity for not life support equipment			
The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.			
Immunity tests	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Product, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF IEC 61000-4-6	3 V rms 150 kHz to 80 MHz outside ISM bands	10 V rms 150 kHz to 80 MHz outside ISM bands	$d = 0,35 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m 80 MHz to 2.5 GHz	$d = 0,35 \sqrt{P}$ 80 MHz to 800 MHz $d = 0,7 \sqrt{P}$ 800 MHz to 2,5 GHz
		w	here P is the maximum output power rating in the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a	Fixed strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To access the electromagnetic environment due to fixed RF transmitters, and electromagnetic site survey should be considered. If the measured field strength in the location in which the Product is used exceeds the applicable RF compliance level above, the Product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Product.		
b	over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.		

Recommended separation distances between portable and mobile RF communications equipment and the not life support equipment			
The Product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Product can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Product as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12 0,	23
0,1	0,38	0,38 0,	73
1	1,2	1,2 2,	3
10	3,8	3,8 7,	3
100	12	12 23	
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the higher frequency range applies.			
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Servicestellen / Servicecenter / Service / Assistenza tecnica / Centro de Servicio

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NOUVAG GmbH
NOUVAG USA Inc.