CGF K1 PC

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 Silfradent S.r.l. Via G. di Vittorio 35/37 47018 S. Sofia (FC) MADE IN ITALY		Manufacturer's data	
	Do not resterilize the device		Keep away from heat sources or direct sun rays
	Caution, consult accompanying documents		The product should be kept dry.
	Do not reuse Disposable device		Consult instructions for use
	LOT AAMM-LLL		REF CGF K1 PC
Batch Code		Device code	
	AAAA-MM		STERILE R
Expiring date		Sterilized using irradiation.	

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Intended use

The CGF K1 PC is a medical device intended to contain the Patient's blood to separate it into its components by centrifugation. It is intended to be used with Phases separator MEDIFUGE.

Fields of use

Dentistry, cosmetic surgery and general surgery. The environmental conditions of use are typical of the dental clinic or surgery.

Device intended exclusively for professional use

Preservation

Recommended temperature for storage: 4–25°C (40–77° F). Take care to maintain the integrity of the sealed pouch.

NOTE:

Avoid direct exposure to the sun. Exceeding the maximum recommended storage temperature may lead to impairment of the tube (e.g. Vacuum loss, etc.).

Lifetime

The device has a lifetime of five years before packaging. The expiring date is indicated both on the data label of the tube and on the sterile pouch.

Disposal

After use the tube must be disposed as a special waste according to legislation in force.

Warning!

Before using the CGF K1 PC always check the integrity of the pouch, if it is damaged do not use it because it is no longer sterile.

The device meets the requirements of CE conformity verified by competent body **0476 Kiwa Cermet**.

CGF K1 PC

Tube intended for separation of blood components for autologous use to be injected / administered to the patient, not for transfusion.

Meaning of the symbols used

In the package of the device and on data label there are these symbols

Silfradent S.r.l. Via G. di Vittorio 35/37 47018 S. Sofia (FC) MADE IN ITALY		Manufacturer's data
	Do not resterilize the device	Keep away from heat sources or direct sun rays
	Caution, consult accompanying documents	The product should be kept dry.
	Do not reuse Disposable device	Consult instructions for use
LOT AAMM-LLL	Batch Code	REF CGF K1 PC Device code
AAAA-MM	Expiring date	STERILE R Sterilized using irradiation.

Instruction for use

Description of the device

The device is a glass test tube of 10 ml capacity vacuum sealed.

Intended use

The CGF K1 PC is a medical device intended to contain the Patient's blood to separate it into its components by centrifugation. It is intended to be used with Phases separator MEDIFUGE.

Fields of use

Dentistry, cosmetic surgery and general surgery. The environmental conditions of use are typical of the dental clinic or surgery.

Device intended exclusively for professional use

Preservation

Recommended temperature for storage: 4–25°C (40–77° F). Take care to maintain the integrity of the sealed pouch.

NOTE:

Avoid direct exposure to the sun. Exceeding the maximum recommended storage temperature may lead to impairment of the tube (e.g. Vacuum loss, etc.).

Lifetime

The device has a lifetime of five years before packaging. The expiring date is indicated both on the data label of the tube and on the sterile pouch.

Disposal

After use the tube must be disposed as a special waste according to legislation in force.

Warning!

Before using the CGF K1 PC always check the integrity of the pouch, if it is damaged do not use it because it is no longer sterile.

The device meets the requirements of CE conformity verified by competent body **0476 Kiwa Cermet**.

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