AS 260

Instructions, Gebrauchsanweisung,
Instrucciones, Instructions, Istruzioni

Contents


Labels and Symbols	4
EN/English.....	8
Notice	8
Overview.....	9
Safety Information	11
Setting Up the Scanner	18
Getting Started	22
Troubleshooting.....	26
Regulatory Information.....	28
Technical Specification	33
DE/Deutsch.....	35
Hinweis	35
Übersicht	36
Informationen zur Sicherheit.....	38
Einrichten des Scanners.....	46
Erste Schritte	50
Fehlersuche	54
Regulatorische Informationen	56
Technische Daten	61
ES/Español	63
Aviso	63
Descripción General	64
Guía de Seguridad	66
Configuración del Escáner	74
Primeros Pasos.....	78
Solución de Problemas.....	82

Información Regulatoria	84
Especificaciones Técnicas	89
FR/Français	91
Avis.....	91
Vue d'Ensemble.....	92
Informations sur la Sécurité	94
Configurer le Scanner	102
Pour Commencer	106
Dépannage	110
Informations Réglementaires	112
Spécifications Techniques	117
IT/Italiano	119
Avviso.....	119
Panoramica	120
Informazioni Sulla Sicurezza	122
Impostazione dello Scanner	130
Per Iniziare.....	134
Risoluzione dei Problemi	138
Informazioni Normative	140
Specifiche Tecniche	145
Contact Information	147

Labels and Symbols

Etiketten und Symbole / Etiquetas y símbolos / Étiquettes et symboles / Etichette e simboli

	EN CE marking applicable for European Union DE CE-Kennzeichnung gilt für die Europäische Union ES Marca CE aplicable para la Unión Europea FR Marquage CE applicable à l'Union européenne IT Marchio CE applicabile all'Unione Europea
	EN Authorized Representative in the European Community DE Bevollmächtigter in der Europäischen Gemeinschaft ES Representante autorizado en la Comunidad Europea FR Mandataire Européen IT Rappresentante autorizzato per la Comunità Europea
	EN Authorized Representative in Switzerland DE Autorisierter Vertreter in der Schweiz ES Representante autorizado en Suiza FR Représentant autorisé en Suisse IT Rappresentante autorizzato in Svizzera
	EN Refer to instruction manual/booklet DE Siehe Gebrauchsinformation ES Consultar el manual de instrucciones/manual FR Se reporter au manuel/notice d'utilisation IT Consultare le istruzioni del Manuale d'Uso
	EN Type BF applied part DE Anwendungsteil vom Typ BF ES Componente aplicado de tipo BF FR Partie appliquée de type BF IT Parte applicata di tipo BF
	EN Manufacturer DE Hersteller ES Fabricante FR Fabricant IT Fabbrikante
	EN Date of Manufacture DE Herstellungsdatum ES Fecha de manufactura FR Date de fabrication IT Data di Produzione

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EN/English

Notice

We recommend that you thoroughly familiarize yourself with this guide to make the most effective use of your system.

The content of this document may change over time. Alliedstar and its subsidiaries are not responsible for any errors in this document or any incidental damages related to its provision, performance, or use. Reproduction of any part of this publication is prohibited without Alliedstar's permission.

This document is originally written in English.

All trademarks and registered trademarks belong to their respective owners.

The AS 260 is intended for professional use only.

Under U.S. federal law, this device can only be sold by a dentist or on a dentist's order.

If any serious incident occurs in relation to the device, the user must report it to Alliedstar and to the competent authority of its Member State in the European Union.

Conventions

The following special notices highlight important information or potential risks to personnel or equipment.

	WARNING: Indicates safety instructions that must be followed precisely to prevent injury to yourself or others.
	CAUTION: Warns of conditions that could cause severe damage or operational issues.
	NOTE: Provides additional information and helpful tips.


Overview

The scanner is designed to acquire 3D models in the followings:

- Upper jaw
- Lower jaw
- Buccal bite registration


Nomenclature

The scanner consists of the following parts.

System Composition

The scanner needs to be connected to a computer. Please note that the computer is not included.

[1] Handpiece
[2] Scanner Tip
[3] Scan Button ● Press to start scanning. Press again to stop scanning.
[4] Mode Switch Button ● Press this button to switch between different modes
[5] Status Indicator ● Rapidly flashing: Scanner is booting. ● Slowly flashing: Scanner is ready but not selected as the active scanner. ● Steady on: Scanner is selected as the active scanner. ● Breathing: Scanner is ready to be connected to the scanning software.
[6] Lens Window
[7] USB Cable
[8] USB Type-C Connector
[9] USB Type-A Connector
[10] Desktop Holder
[11] Wall Mount Holder
[12] Software

Holder

Place the scanner in the holder when you are not using it.

NOTE: The scanner goes into sleep mode when put into the holder or left idle for 10 seconds. To use it again, take it out of the holder or pick it up and press any button on the scanner.

Safety Information

Intended Use

The AS 260 is a digital optical scanning device used to obtain digital impressions of hard and soft tissues such as teeth, gingival, and mucous membrane by means of oral scanning, for oral restoration and orthodontic treatment of malocclusion.

The AS 260 could be used for both adults and children in clinical practice.

Clinical Benefits and Performance Characteristics

The AS 260 benefits a dental practice by enabling practitioners to acquire digital impressions with the quality and accuracy required for digital CAD/CAM dental applications. The actual performance of the device is dependent on the user's training and operating execution. The user is solely responsible for the accuracy, completeness, and adequacy of the acquired data.

Contraindications

None.

Warnings and Safety Instructions

DANGER OF ELECTRIC SHOCK: This is an electrical unit. **DO NOT** expose it to water spray. Such action can cause an electric shock or a malfunction of the unit.

CAUTION: All known residual risks, contraindications, or undesirable side effects are listed in this guide. If any serious incident occurs in relation to the device, you must report it to Alliedstar and to the competent authority of your Member State in the European Union.

WARNING

AS 260

- Before using the scanner, you must read and understand all safety information.
- Use the scanner only within hospitals and professional healthcare facilities. Do not operate it near high-frequency surgical equipment or inside the RF-shielded room of an ME system for magnetic resonance imaging, where electromagnetic disturbances are intense.

- Before operating the scanner, inspect the scanner and its accessories for any rough surfaces, sharp edges, or protrusions that could pose a safety hazard.
- You are responsible for the scanner's operation and maintenance. Training is required to use the scanner safely.
- Do not place objects within the scanner's operational field.
- Ensure the scanner is turned off when not in use.
- Do not use the scanner in oxygen-rich environments. It is not intended for use with flammable anesthetics or agents.
- Avoid pulling or twisting the cable.
- Do not drop the scanner or its accessories.
- Do not attempt to sterilize the scanner.
- Keep the scanner away from water sprays, and do not immerse it in water or disinfectants.
- Avoid subjecting the scanner to heavy vibrations.
- Avoid exposing the scanner to ultraviolet radiation directly. The scanner is not designed for ultraviolet disinfection.
- Avoid staring directly at the LED emission window.
- Avoid removing the cover from any part of the scanner. The device has no user-serviceable components. If repairs are needed, reach out to a qualified Alliedstar service technician.
- Do not use any cables other than those provided with the scanner. Using different cables may damage the scanner and compromise its safety features and EMC performance.
- Keep equipment that does not comply with IEC-60601 standards at least 1.5 meters from the patient.
- If the equipment malfunctions, switch off, post an "Out of Service" notice, and contact a certified Alliedstar technician.
- The use of components, accessories, cables, or spare parts not specified or supplied by the manufacturer of this equipment can compromise the scanner's safety features. This may lead to higher electromagnetic emissions, reduced electromagnetic immunity, and improper equipment functioning.
- Do not modify this equipment in any way.
- Do not connect additional power strips or extension cords to the system.
- The applied part may reach temperatures up to 43 °C (109.4 °F). To prevent overheating, avoid prolonged use.
- To completely disconnect the device from the main power supply, unplug the USB connector from the USB port.
- DO NOT maintain or service this equipment while it is in use with the patient.

- Connection of the PEMS (Programmable Electrical Medical System) to an IT NETWORK that includes other equipment could result in risks to patients, operators, or third parties. The responsible organization should identify, analyze, evaluate, and control these risks.
- Patients with oral mucosal disease, mental illness, severe respiratory disease, asthma, Parkinson's disease, hyperactivity disease are forbidden.
- Patients with moderate or severe opening limitation should use it with caution.

Computer

- DO NOT place any equipment which does not comply with IEC 60601-1 in the immediate vicinity of the patient. Leave at least 1.5 meters distance between the patient and the equipment.
- The scanner is only intended to be connected to a computer that is at least IEC 60950 / IEC 62368, or equivalent standards certified. Connecting the scanner to other equipment may be hazardous.
- Refer to your computer's installation guide for details on the data processing system, computer, and monitor. Ensure there is enough space around the computer to allow for proper ventilation.
- Place the screen so that light reflections from both internal and external sources are minimized to achieve the best image quality and ensure comfortable viewing.

Disposal

This device contains materials and chemicals commonly used in electrical and electronic products. Improper disposal at the end of its life can lead to environmental pollution. Therefore, do not discard it with regular household waste. Please take it to a designated electronic waste disposal or recycling facility. For more information on how to dispose of electronic waste, please get in touch with your local authorities. Dispose of the Scanner Tips according to standard operating procedures or local regulations for the disposal of contaminated medical waste. For additional Scanner Tips, contact your dealer.

Cleaning, Disinfecting, Sterilizing

Regularly perform the maintenance tasks outlined below for your scanner and accessories.

Carefully follow the instructions when preparing the scanner for use to ensure the highest level of hygienic safety for patients.

To ensure maximum hygienic safety for the patient and to minimize the risk of cross-contamination, carefully perform the following maintenance activities on your scanner and accessories.

After each patient:

- Clean and disinfect the scanner. (See 'Clean and Disinfect the Scanner')
- Clean the Scanner Tip, and then perform Autoclave Sterilization (See 'Clean the Scanner Tip' and 'Sterilize the Scanner Tip').

There are 3 product models of Scanner Tip:

Product Model	Size	UDI-DI	Manual Cleaning	134°C Sterilization
TP202	Large	(01)06973993441034	Yes	Yes
TP203	Small	(01)06973993441041	Yes	Yes
TP204	Large	(01)06973993441386	Yes	Yes

Clean and Disinfect the Scanner

General Warnings

WARNING

- When handling the scanner, read and follow all warnings and personal protection guidelines in the sanitizer's Safety Data Sheet (SDS).
- Gloves must be worn when cleaning and disinfecting the scanner.
- After each patient, disinfect the scanner using a recommended medium-level disinfectant wipe with tuberculocidal properties.
- Do not use disinfectants containing phenol or iodine-based substances, as they may damage the scanner's surface coating.
- Do not place the scanner in a pressure steam sterilizer or immerse it in water or disinfectant solutions.
- Excessive liquid can harm the scanner.
- Avoid using cotton balls, cloth pieces, or paper towels soaked in disinfectant to sanitize the scanner.

Clean the Scanner

If visible blood or body fluid is on the scanner, clean it before proceeding with disinfection.

To clean the scanner, follow the below steps:

- Dampen a lint-free cloth with warm water (do not soak).
- Use the dampened lint-free cloth to wipe away blood or body fluids.

Disinfect the Scanner

The scanner must be thoroughly sanitized after use for each patient.

To properly disinfect the scanner, make sure to follow the disinfectant manufacturer's guidelines regarding the required contact time.

CAUTION: The scanner must be cleaned before disinfection if it is visibly soiled. Refer to "Cleaning the Scanner."

To disinfect the scanner, proceed with the following steps:

- Detach the Scanner Tip.
- Clear any visible debris (refer to "Cleaning the scanner").
- Use an intermediate-level disinfectant wipe that is commercially available. Be sure to follow the manufacturer's instructions for the necessary contact time.

Recommended disinfectant wipes: CaviWipes

WARNING: Using an unapproved disinfectant may damage the scanner.

- Thoroughly wipe down all scanner surfaces, ensuring no liquid enters any gaps, air outlets, or pinholes.

WARNING: Do not rinse the scanner.

- Let it air dry thoroughly.
- Once dry, use a clean, lint-free cloth dampened with water to remove any remaining disinfectant from the scanner's surface.

Clean and Sterilize the Scanner Tip

WARNING

- Always wear gloves when handling a contaminated Scanner Tip.
- Review and follow the safety and personal protection instructions in the detergent manufacturer's Safety Data

Sheet (SDS).

- Avoid soaking the Scanner Tip in disinfectant for extended periods.
- Ensure the scanner tip is completely dry before mounting it to the scanner.
- Do not use an ultrasonic cleaner to clean the Scanner Tip.
- Do not soak the Scanner Tip with alcohol-based disinfectants.

Do not soak the Scanner Tip with alcohol-based disinfectants. Clean the Scanner Tip

To manually Clean the scanner Tip, follow these steps:

- Rinse the scanner tip under running water to remove excess debris (approximately 2 minutes).
- brush all surfaces with an enzymatic detergent solution, like Metrex EmPower.
- Rinse thoroughly under fresh, running water for a full 2 minutes.
- Inspect the scanner tip. If it remains dirty, go through the steps above again.
- Clean the mirror inside the Scanner Tip using a lint-free cloth or lens tissue to remove dust.
- Dry the scanner tip's mirror carefully with a lens tissue or lint-free cloth.

Sterilize the Scanner Tip

Scanner Tips shipped non-sterile. You must sterilize them before use.

CAUTION: For TP202 and TP203, if you limit the exposure time at 134°C to not more than 6 minutes, you can autoclave the Scanner Tip up to 60 cycles.

CAUTION: For TP204, if you limit the exposure time at 134°C to not more than 6 minutes, you can autoclave the Scanner Tip up to 180 cycles.

To sterilize the Scanner Tip, follow these steps:

- Put the Scanner Tip into a sealed steam sterilization pouch.
- Place the Scanner Tip in a steam autoclave for sterilization.
 - ◆ Exposure temperature should be set as 134°C.
 - ◆ Exposure time should exceed 3 minutes.
 - ◆ Exposure time should not exceed 6 minutes.

Precautions before Use

Perform the following activities on your product and accessories before use.

Visually Inspecting the Scanner for Damage

Visually inspect the scanner for damage or signs of deterioration by doing the following:

- Inspect the scanner's Lens Window.
- Inspect around the scanner buttons and contact points.

If damage is noted, do not use the scanner, and contact your representative.

Visually Inspecting the Scanner Tips

Visually inspect the Scanner Tips for signs of deterioration by doing the following:

- Verify that the Scanner Tip is not damaged, and its components are not detached.
- Verify that the Scanner Tip mirror does not have any smudges or scratches on it.
- If deterioration is noted, replace the Scanner Tip.

WARNING

The Lens Window on the scanner is a delicate optical component. Mount the lens window protective cover to protect the Lens Window from damage and dirt when the scanner is not in use.


The mirror in the Scanner Tip is a delicate optical component. Its clean and undamaged surface is critical to scan quality.

If you see poor scan quality or an unclear video preview in the software, clean the mirror of the Scanner Tip and the Lens Window using a microfiber cleaning swab, applying ethanol that is free of impurities.


Setting Up the Scanner

To set up the scanner, follow these steps:


- Go to <https://support.allied-star.com/download> and download the installation file and instructions according to the product model.
- Double-click the software installation file.
- Choose a language from the drop-down list and click Ok to install.
- Follow the instructions on the screen to complete the installation.
- Insert the USB Type-C Connector. Then push up the cable thumb push sleeve to secure the USB Type-C Connector.

- To disconnect the USB Type-C Connector, use thumb and index finger to pinch the area with red circle and pull the cable thumb push sleeve off. Then detach the USB Type-C Connector from its socket.

- Remove the lens window protective cover on the Lens Window.
- Firmly mount the Scanner Tip onto the Lens Window.

- Insert the Type-A Connector to a USB 3.0 Type-A port on the computer.

CAUTION: Make sure the scanner is connected to the USB 3.0 port. If connected to a USB 2.0 port, the scanner may not work properly.

- The scanner will automatically turn on, the status indicator will flash rapidly, then the Status Indicator will stay on once it's connected to the computer and ready to scan.


NOTE: When the scanner is connected to the software, the indicator stays in solid on.

Using the Holder

Installing the Desktop Holder

To use the holder as a Desktop Holder, follow these steps:


- Select a clean surface area.
- Place the scanner in the Desktop Holder, choose an appropriate length and fix the USB cable to cable clips.


Installing the Wall Mount Holder

To use the holder as a Wall Mount Holder, follow these steps:


- Select a flat and secure area with easy access.
- Insert screws through the holes in the holder rail and wall mount pad in turn to affix it to a solid surface.

- Slide the Wall Mount Holder into the holder rail.

- Slide the scanner into the Wall Mount Holder

NOTE: If the Wall Mount Holder is not properly installed, there is a risk that the holder can fall off the wall, resulting in damage to the scanner.

Getting Started

Accessing the Software User Interface

To access the software user interface, follow these steps:

- Double-click the software icon on the desktop.
- The software will automatically pop up the login window.
- If you don't have an account, click **Create an account** to register the organization, and complete the email verification.
- Type your account information in the username and password fields, click the **Log in** button.
- The **Cases** page will be shown up.
- Click the **New Case** button, then input patient name and click the **OK** button.
- Click the **Scan** button.
- If the scanner is not activated, the device activation dialog will be displayed. Follow the instructions on the screen to complete device activation.
- Click the Option menu button  and select **Preferences**.
- Customize the configuration options.
- You can now start acquiring 3D models.

Audio Prompts

Scanning Sound

When the **Enable scanning sound** option is enabled, your computer will play a continuous sound when you are successfully scanning. If the sound stops, it means the scan has stopped. If you need to continue, please return to the previous scanning area until the scanner resumes scanning and your computer plays a continuous sound. When a buccal bite registration is successfully scanned, your computer will also play a short sound.

Warning Sound

When the **Enable warning Sound** option is enabled, if the cumulative scanning time of the current case exceeds the recommended threshold and your computer may not be able to maintain the peak scanning performance, your computer will play a short warning sound.

NOTE: If your computer does not have speakers, these options will not take effect.

Preparing the Teeth

- If there is a preparation area, retract the gingiva by gingival restriction cords. And extract the cords just before scanning the preparation.
- Before starting the scan, dry the teeth thoroughly.
- During the scan, re-dry the teeth moderately.

Preparing the Scanner


The Scanner Tip attached to the scanner provides sanitary shield for patients. Always disinfect the body of the scanner and perform sterilization on the Scanner Tip after each use.

CAUTION: Scanner Tip received from the manufacturer are NOT sterilized. You must sterilize them before the first use.


CAUTION: Avoid any liquid from leaking into the air outlet near the Lens Window or the air inlet at the rear of the scanner (see the figure below), otherwise the scanner may be damaged.

To prepare the scanner, follow these steps:

- Make sure the Lens Window at the base of the scanner is clean by wiping it with a moist, lint-free cloth or lens tissue.
- Slide the Scanner Tip onto the scanner as shown below.


Starting Scanning

To start scanning, place the Scanner Tip on the surface of the tooth to stabilize the scanner and press the Scan button. Wait until a 3D model appears on the 3D model display screen, and then slowly move it along the arch at 0-5mm from the teeth.

Scanning Approach

The recommended scanning method is to start with a molar, since it has greater details for easier identification. Change the scanning angle to less than 60 degrees during scanning to allow the surfaces to overlap, if the overlap is too small, the alignment may be lost.


Scanning Protocol

The recommended scanning protocol consists of 3 sweeps: occlusal, lingual and buccal to ensure good data coverage of all surfaces.

It is recommended to start the first sweep from the occlusal surface, you should start with the first molar.

The second sweep can scan both the lingual and buccal sides, and the third scan covers the opposite side of the

second sweep.

Troubleshooting

AS 260 Troubleshooting Instructions

Problem Description	Action
No connection between the scanner and the computer.	Detach the USB Type-C Connector from the Handpiece, and insert again. If the problem is not solved, replace the cable.
The 3D model shows discrepancies and overlapping areas.	Perform the scan using the suggested scanning method. Use the Cut tool to eliminate mismatched data and excess tissue, then conduct a rescan.
Following the buccal bite registration, a gap or an overlap may appear between the upper and lower jaws.	● Adjust the Occlusal pressure level in Preferences, then disable and enable Occlusal pressure adjustment. ● Delete the incorrect bite view and rescan.
A decline in precision or poorly stitched images may be noticed during acquisition.	Clean the lens window at the scanner's base using a moist, lint-free cloth or lens tissue to remove dust or water stains. Ensure the Scanner Tip is securely attached, and check that the preview video does not have dark edges or other dirt.
Reconstructing metallic preparations can sometimes be challenging.	Adjust the scanner's position—such as its distance or angle—and capture more of the area during the scan. Move the surgical light away from the patient to decrease light scatter. Turn on Shining Surface.
The Scanner Tip is installed but not detected. No live video is displayed, and the Scanner Tip not detected icon is displayed at the bottom-right of the interface.	After removing the scanner tip, wipe the Scanner Tip Detection Sensor to ensure it is free from dirt or debris. Reattach the Scanner Tip and verify that it is securely and firmly connected to the scanner.
The inner surface of the scanner's lens window at its base shows signs of fogging.	Attach a fully dry Scanner Tip to the scanner, place the scanner in the Holder or desk, and wait for the fogging to clear. If the fogging persists after 24 hours, contact your local service provider for help. Ensure the Scanner Tip is completely dry before attaching it, and avoid using disinfectant-soaked cloths to clean the scanner.
The scanner does not cast light, and a static preview image is displayed in the video preview screen.	Ensure that your computer meets the Minimum System Requirements. Close the software, and then launch the software again.
After stopping scanning, the scanner is disconnected and cannot be	Try the following steps in order. Check if the scanner is turned off. Press any button on the scanner or

automatically reconnected.	take it out of the holder to turn on the scanner again. The scanner startup process may take several seconds. Disconnect the USB connector of the scanner from the computer, and then reconnect it to the computer. The scanner will automatically power on and connect.
----------------------------	--

Regulatory Information

The AS 260 complies with the following regulations:

- (EU) 2017/745 Medical Device Regulation (MDR).
- FDA Center for Devices & Radiological Health CDRH - Title 21 CFR 872.3661 (USA).
- Medical Devices Regulations (Canada).
- Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment, 2011/65/EU Annex II and its amendment Directive (EU) 2015/863.

Compliance with European and International Standards

EN 60601-1 / IEC 60601-1: Medical Electrical Equipment, Part 1: General Requirements for Basic Safety and Essential Performance

ANSI / AAMI ES 60601-1: Medical Electrical Equipment - Part 1: General requirements for basic safety and essential performance

CAN / CSA-C22.2 No. 60601-1: Medical Electrical Equipment - Part 1: General requirements for basic safety and essential performance

EN 60601-1-2 / IEC 60601-1-2: Medical Electrical Equipment, Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – Requirements and tests

EN 80601-2-60 / IEC 80601-2-60: Medical electrical equipment — Part 2-60: Particular requirements for the basic safety and essential performance of dental equipment

EN 62471 / IEC 62471: Photobiological safety of lamps and lamp systems: Equipment classification, requirements, and User's Guide

EN ISO 17664: Processing of health care products – Information to be provided by the medical device manufacturer for the processing of medical devices

EN ISO 17665-1: Sterilization of health care products — Moist heat — Part 1: Requirements for the development, validation, and routine control of a sterilization process for medical devices

EN 60601-1-6 / IEC 60601-1-6: Medical Electrical Equipment, Part 1-6: General requirements for basic safety and essential performance - Collateral Standard: Usability

EN 62366-1/ IEC 62366-1: Medical devices - Application of usability engineering to medical devices

EN 62304 / IEC 62304: Medical device software - Software life cycle Processes

EN ISO 10993: Biological evaluation of medical devices

EN ISO 14971: Medical devices - Application of risk management to medical devices

EN ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied - Part 1: General requirements

EN ISO 20417: Medical devices — Information to be supplied by the manufacturer

ISO 9687: Dentistry - Graphical symbols for dental equipment

AAMI TIR 12: Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers

AAMI TIR 30: A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices

Classification in Accordance with EN / IEC 60601-1

Type of protection against electric shock: N/A

Degree of protection against electric shock: Type BF Applied Part

Mode of operation: Continuous operation

Flammable anesthetics: Not suitable for use in the presence of flammable anesthetics or a mixture of flammable anesthetics with air or oxygen or nitrous oxide.

Conformity with EN/ IEC 60601-1-2

IEC 60601-1-2 EMC requirements and tests, Medical Electrical Equipment including CISPR 11 Group 1, Class B.

Electromagnetic Compatibility Precautions

Medical electrical devices need specific measures to ensure electromagnetic compatibility (EMC).

Medical devices should be set up and used following the EMC guidelines detailed in this document.

Even equipment that meets CISPR emission standards may disrupt communication with the AS 260.

WARNING: To maintain the AS 260's optimal performance, ensure that all portable and mobile RF communication devices—including peripherals like antenna cables and external antennas—are used at a distance of no less than 30 cm (12 inches) from any part of the system, including the cables specified by the manufacturer. Using such equipment closer than this recommended distance may result in decreased performance of the AS 260.

Guidance and Manufacturer's Declarations

Guidance and Manufacturer's Declaration - Electromagnetic Emissions

The AS 260 is intended for use in the electromagnetic environment specified below. The customer or user of the AS 260 should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1 Class B	The AS 260 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity for Equipment and Systems

The AS 260 should be operated only within the electromagnetic environments described below for optimal performance. Customers or users must ensure that the device is used under these specified conditions.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	For optimal conditions, floors should be constructed from wood, concrete, or ceramic tile. If synthetic materials are used to cover the floors, please ensure that the relative humidity is maintained at a minimum of 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	For optimal conditions, power-frequency magnetic fields should be maintained at levels commonly found in standard commercial or hospital environments.
Radiated RF IEC 61000-4-3	3 V/m 80MHz – 2.7GHz	3 V/m 80MHz – 2.7GHz	In a professional healthcare setting: Warning: To prevent potential interference, portable RF communications devices (including accessories like antenna cables and external antennas) must be kept at least 30 cm (12 inches) away from any part of the AS 260, including any cables specified

			by the manufacturer. Failure to maintain this distance could reduce the equipment's performance.
--	--	--	--

The AS 260 meets the immunity requirements for proximity fields from RF wireless communication devices, in line with the test levels outlined below according to the IEC 60601-1-2 standard. It is the customer or user's responsibility to ensure that the AS 260 is operated in a suitable environment.

Test Frequency (MHz)	Band (MHz)	Immunity Test Levels
385	380 – 390	Pulse modulation 18Hz, 27V/m
450	430 – 470	FM, ± 5 kHz deviation, 1 kHz sine, 28V/m
710	704 – 787	Pulse modulation 217Hz, 9V/m
745		
780		
810	800 – 960	Pulse modulation 18Hz, 28V/m
870		
930		
1720	1700 – 1990	Pulse modulation 217Hz, 28V/m
1845		
1970		
2450	2400 – 2570	Pulse modulation 217Hz, 28V/m
5240	5100 – 5800	Pulse modulation 217Hz, 9V/m
5500		
5785		

Accessories

Using cables or accessories other than those specified—unless they are manufacturer-approved replacements for internal components—could lead to increased emissions or reduced immunity of the medical device.

Other Equipment

WARNING: Avoid placing this equipment next to or stacking it with other devices, as this may cause improper functioning. If such placement is necessary, monitor this equipment and the other devices to ensure they operate correctly.

Technical Specification

Product Model

AS 260

AS 260 Technical Specifications

Components	Technical Specifications
Weight	Handpiece: 175g Tip Large - Sterilization (TP202): 14g Tip Small - Sterilization (TP203): 12g Tip Large - Sterilization (TP204): 14g
Color	3D Full Color
Connectivity	USB 3.0
Power Source	USB 3.0 5V, 900mA
Field of View	Tip Large - Sterilization (TP202): 16mm x 14mm Tip Small - Sterilization (TP203): 12mm x 12mm Tip Large - Sterilization (TP204): 16mm x 14mm
Minimum System Requirements of Laptop	Processor: 10th Generation Intel® Core™ i7-10750H or AMD Ryzen™ 7 4700U Memory: 16 GB Graphics Card: NVIDIA GeForce GTX 1650 4GB Disk: 512 GB SSD Operating System: Win 10 (build 18362+) / Win 11, 64 bit Others: USB 3.0 Port Optional: Touch Screen
Minimum System Requirements of Desktop	Processor: 9th Generation Intel® Core™ i5-9600 or AMD Ryzen™ 5 2600 Memory: 16 GB Graphics Card: NVIDIA GeForce GTX 1650 4GB Disk: 512 GB SSD Operating System: Win 10 (build 18362+) / Win 11, 64 bit Others: USB 3.0 Port Optional: Touch Screen

CAUTION: It is MANDATORY to check that your system configuration is compatible with the computer system requirements for the AS 260 software.

Environmental Requirements

Components	Environmental Requirements
Operating Temperature	15°C ~ 30°C
Transport and Storage Temperature	-10°C ~ 50°C
Operating Relative Humidity	10% ~ 65% RH
Transportation and Storage Relative Humidity	10% ~ 95% RH
Operating Atmospheric Pressure	70 ~ 106 KPa
Transportation and Storage Atmospheric Pressure	60 ~ 106 KPa

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